



ADVANCE CARE **IP** Air Bubble Call Pendant

The Advance Care Air Bubble Call Pendant has been designed with the needs of 'restricted movement' residents in mind - for example those with arthritis. Only the slightest pressure is required from the resident to activate the device, thanks to the cleverly designed pneumatic activation. The device is easy to hold as it is very light and has a clip for easy attachment to resident's bed or clothing to keep them within easy reach.

As the Air Bubble Call Pendant is plugged into the Nurse Call Button via the tough, fully moulded 6.35mm jack into the pendant socket on the front of the plate, the Call Pendant essentially acts as an extension of the Nurse Call Button, allowing for the easy access of a Call Point without having to strain reaching the wall.

Although an extension of the Nurse Call Button, when activated, the Air Bubble Call Pendant has its own easily identifiable alert, such as 'Pendant Room 1'.

When a call is made, a high visibility reassurance light illuminates the 'CANCEL' button on the Call Button plate and the Call Button may also 'beep' periodically (if activated), providing reassurance to the user.

The Nurse Call alert is also immediately displayed at the Nurse Station and sent to associated paging devices.

The Advance Care IP Nurse Call System includes a 'Pendant Removed' alarm feature which generates an alert if the pendant is accidentally pulled out or incorrectly plugged in, a beneficial feature with resident safety as key.

The pendant itself is totally waterproof, completely safe in oxygen rich environments and is UL1069 listed.

Available with an 8 or 12 foot cable, the Air Bubble Call Pendant integrates with most Nurse Call Systems and is a valued addition to any Advance Care IP Nurse Call installation.

- ### Features
- Soft pneumatic activation
 - Suitable for residents with restricted movement
 - Integrates with most Nurse Call Systems
 - 6.35mm industry standard Pendant Plug
 - Safe in Oxygen rich environments; UL1069 rated
 - Self-configured alarm notification
 - Bed clip as standard
 - Cord length options



Code	Description
AC-PA/8	Air Bubble Call Pendant (2400mm/8 foot cord length) - white or grey
AC-PA/12	Air Bubble Call Pendant (3700mm/12 foot cord length) - white or grey

Tender Specifications

Supply and install minimum dexterity Call Pendants where indicated on the drawings.

The Pendant will be ergonomically designed for ease of grip by residents to activate calls.

The Pendant is to be waterproof and completely safe in oxygen rich environments (UL1069 listed).

The Pendant shall be tamper proof against disassembly by residents, children, staff and intentional or unintentional impact. The supplier shall provide a two year warranty against damage to the pendant body.

The Pendant Call Button will be:

- Activated through a momentary press of the Air Bubble
- Easy to use for residents with restricted movement and require the slightest press to activate
- Provide the resident with a positive tactile response when pressed

The Pendant's cord will be supplied with an integrated 8 or 12 foot multi-stranded cable (as specified in the schedule of prices), double insulated white PVC cable with an injection moulded right angle 6.35mm phono plug.

The cord will have a large, non-interfering metal bed linen clip attached to the pendant cable. The clip shall be removable and replaceable.

Removal of the Pendant, either intentionally or by accident, will result in a call being activated stating 'Pendant Removed'. The pendant removal alarm can only be reset by reinserting a Pendant Plug or by inserting Dummy Plug.

Design Specifications

Feature	Description
Height	25mm
Width	40mm
Depth	38mm
Weight	190 grams
Connectors	6.35mm Stereo Phono Plug
Temperature	-20°C to 50°C
Humidity Range	0 to 90%
Rating	UL1069

Compliance

Designed to comply with:

- AS/NZS 3811 Hard Wired Patient Alarm Systems
- AS/NZS 2999 Alarm Systems for the Elderly and Other Persons at Risk
- AS/NZS 3000 Wiring Rules
- AS/NZS CISPR22 (Class B) Electromagnetic Compatibility (EMC)

© 2015 Digital Telecommunication Systems (WA) Pty Ltd

All rights are reserved.

Digital Telecommunication Systems (WA) Pty Ltd reserves the right to make changes in specifications and/or discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Authorised Distributor

